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# Table of contents

---

1. Scale models.....	28
2. Scope of delivery.....	28
3. Warnings.....	29
4. General notes .....	30
4.1 Technical data.....	30
4.2 Intended use .....	32
4.3 Classification.....	32
4.4 Safety instructions.....	33
4.5 Cleaning.....	34
4.6 Maintenance and service .....	34
4.7 Guarantee / Warranty / Liability.....	34
4.8 Battery and accumulator disposal .....	35
4.9 Disposal of the scale.....	35
5. Labelling.....	36
6. Assembly and installation - chair scales.....	37
6.1 Unpacking the scale.....	37
6.2 Mounting the scale .....	38
7. Assembly and installation - wheelchair scales.....	39
7.1 Mounting the scale.....	39
7.2 Setting up the scale .....	40
7.3 Transport of the scale.....	40
8. Assembly and installation - person scales.....	41
8.1 Setting up the scale .....	41
9. Product description.....	42
9.1 Display .....	42
9.2 Operating elements .....	43
10. Basic functions .....	44
10.1 Switching on the scale.....	44
10.2 Weighing.....	44
10.3 Tare.....	44
10.4 Zeros .....	44
10.5 BMI function .....	45
10.6 Switching off the scale .....	45
11. Communication protocols.....	45
12. Malfunctions - causes and elimination .....	46
13. Note on electromagnetic compatibility .....	47

# 1. Scale models

## Conformity assessed

## Not conformity assessed

These instructions for use apply to the following scale models:

- ▶ 6702.01.402 Chair scale
- ▶ 6708.01.502 Wheelchair scale
- ▶ 6710.01.002 Personal scale

- ▶ 6802.01.402 Chair scale
- ▶ 6808.01.502 Wheelchair scale
- ▶ 6810.01.002 Personal scale

# 2. Scope of delivery

Chair scale 6702:

Article number	Item number
Chair weigher	6702.01.402
Power supply unit	618.020.102
Instruction manual	470.051.221

Chair scale 6802:

Article number	Item number
Chair weigher	6802.01.402
Power supply unit	618.020.102
Gebrauchsanweisung	470.051.221

Wheelchair scales 6708:

Article number	Item number
Wheelchair scale	6708.01.502
Power supply unit	618.020.102
Instruction manual	470.051.221

Wheelchair scales 6808:

Article number	Item number
Wheelchair scale	6808.01.502
Power supply unit	618.020.102
Instruction manual	470.051.221

Personal scale 6710:

Article number	Item number
Personal scales	6708.01.502
Power supply unit	618.020.102
Instruction manual	470.051.221

Personal scale 6810:

Article number	Item number
Personal scales	6810.01.002
Power supply unit	618.020.102
Instruction manual	470.051.221

### 3. Warnings

---

- ▶ In the event of a malfunction, contact the dealer or the manufacturer. Unauthorised modifications or repairs may damage your scale and void the manufacturer's warranty.
- ▶ The indicator and scale must not get wet. Liquids (e.g. water) can cause damage. Use a dry cloth, e.g. a towel, to dry the indicator.
- ▶ Disconnect the power supply to this unit before any installation, cleaning or maintenance. Failure to do so could damage the unit.
- ▶ If the scale is not to be used for a longer period of time, please remove the power supply.
- ▶ Avoid stacking materials on the display or placing weights on the display. This can cause damage.
- ▶ Place the scale on a firm, stable and level surface so that accurate measurement results can be guaranteed. If the surface is soft or inclined, the measurement results will not be representative.
- ▶ Do not connect the indicator to unstable power sources.
- ▶ Use only original equipment. Use of other makes may cause damage to the scale.
- ▶ Do not touch the interface and the patient at the same time.

Bump, crush, fall or trip hazard

Follow the instructions for use

## 4. General information

### 4.1 Technical data

	Chair scale 6702	Chair scale 6802
Scale type	Single range scale	Single range scale
Maximum load	300 kg	300 kg
Minimum load	2 kg	2 kg
Digit increment	100 g	100 g
Dimensions (W x D x H)	850 x 650 x 970 mm	850 x 650 x 970 mm
Product weight	18 kg	18 kg
Power supply	Battery or mains operation Adapter 120VAC-9VDC-50Hz / 230VAC-9VDC-50Hz Use only power supplies that are approved for medical applications applications and comply with the EN 60601 standard.	Battery or mains operation Adapter 120VAC-9VDC-50Hz / 230VAC-9VDC-50Hz
Calibration class	Calibration class III, MPG class 1 m	
Working temperature	+ 10 °C bis + 35 °C	+ 10 °C bis + 35 °C
Storage temperature	0 °C bis + 50 °C	0 °C bis + 50 °C
Humidity	85 % (non-condensing)	85 % (non-condensing)
Air pressure	Use at atmospheric pressure	Use at atmospheric pressure

	Wheelchair scales 6708	Wheelchair scales 6808
Scale type	Single range scale	Single range scale
Maximum load	300 kg	300 kg
Minimum load	2 kg	2 kg
Digit increment	100 g	100 g
Dimensions (W x D x H)	1.250 x 1.000 x 1.150 mm	1.250 x 1.000 x 1.150 mm
Product weight	46 kg	46 kg
Power supply	Battery or mains operation Adapter 120VAC-9VDC-50Hz / 230VAC-9VDC-50Hz Use only power supplies that are approved for medical applications applications and comply with the EN 60601 standard.	Battery or mains operation Adapter 120VAC-9VDC-50Hz / 230VAC-9VDC-50Hz
Calibration class	Calibration class, MPG Klasse 1 m	
Working temperature	+ 10 °C bis + 35 °C	+ 10 °C bis + 35 °C
Storage temperature	0 °C bis + 50 °C	0 °C bis + 50 °C
Humidity	85 % (non-condensing)	85 % (non-condensing)
Air pressure	Use at atmospheric pressure	Use at atmospheric pressure

## 4. General information

	Personal scale 6710	Personal scale 6810
Scale type	Single range scale	Single range scale
Maximum load	250 kg	250 kg
Minimum load	2 kg	2 kg
Digit increment	100 g	100 g
Dimensions (W x D x H)	375 x 375 x 78 mm	375 x 375 x 78 mm
Product weight	10,5 kg	10,5 kg
Power supply	Battery or mains operation Adapter 120VAC-9VDC-50Hz / 230VAC-9VDC-50Hz Use only power supplies that are approved for medical applications applications and comply with the EN 60601 standard.	Battery or mains operation Adapter 120VAC-9VDC-50Hz / 230VAC-9VDC-50Hz
Calibration class	Calibration class III, MPG class 1 m	
Working temperature	+ 10 °C bis + 35 °C	+ 10 °C bis + 35 °C
Storage temperature	0 °C bis + 50 °C	0 °C bis + 50 °C
Humidity	85 % (non-condensing)	85 % (non-condensing)
Air pressure	Use at atmospheric pressure	Use at atmospheric pressure

## 4. General information

---

### 4.2 Intended use

These scales are used to determine the weight of persons and may only be used for this purpose. Any other use is prohibited. These scales are intended for use in legal metrology. The scales comply with the provisions of the applicable European directives and regulations.

The manufacturer's medical product advisor must be notified in writing of any faults that could endanger the patient or lead to the falsification of measurement results. These instructions for use are an integral part of the devices. Precise observance of these instructions is a prerequisite for the intended use and correct operation of the devices.

Operation is only to be carried out by qualified personnel. Operation by the patient is not intended.

Please ensure that the scales are not subjected to permanent loads. In addition, shocks and overloading of the scales in excess of the specified maximum load must be avoided at all costs. This can damage the scales

### 4.3 Classification

Electrical protection class II (protective insulation, no protective ear connection).

#### Application part (6702 & 6708 & 6710)

The scale is a class 1 medical device with a measuring function.

The medical device is a type B applied part and is intended for direct contact by the patient. The leakage currents correspond to the classification of type B applied parts..

## 4. General information

---

### 4.4 Safety instructions

These scales are designed for medical use. The user must be familiar with the operation of the scale. Please read the information in the instruction manual carefully before using the scale. It contains important instructions for the installation, the intended use and the maintenance of the device. The manufacturer is not liable if the following instructions are not observed. Improper installation will invalidate the warranty.

- ▶ The scale is not explosion-proof and must not be used in critical or potentially explosive environments.
- ▶ When using electrical components under increased safety requirements, the corresponding regulations must be observed.
- ▶ Electrical connection conditions must comply with the values printed on the power supply unit
- ▶ The scale is designed for indoor operation. Observe the permissible ambient conditions for operation (see technical data). The scale meets the requirements for electromagnetic compatibility. Exceeding the maximum values specified in the standards must be avoided
- ▶ The scale must be installed in such a way that it can be disconnected from the mains easily and at any time. The accessibility of the socket should be guaranteed at all times.
- ▶ When moving the scale, the scale must be raised. The scale must never be pushed. Then check the levelling and adjust it to the new situation if necessary.
- ▶ The scale is a measuring instrument. Draughts, vibrations, rapid temperature changes and sunlight can affect the weighing result. High humidity, vapours, aggressive liquids and heavy soiling must be avoided.
- ▶ This unit is interference suppressed in accordance with the applicable EC directive. However, under extreme electrostatic and electromagnetic influences, e.g. when operating a radio or mobile phone in the immediate vicinity of the unit, the display value may be affected. After the interference has ceased, the product can be used again as intended; if necessary, it must be switched on again.
- ▶ Armrests and footrests of the chair scale 6702 are not standing aids. They may only be used as support surfaces for arms and feet.
- ▶ The chair scales 6702 and 6802 are not a transport aid for persons or objects. There is a risk of injury.

## 4. General information

---

### 4.5 Cleaning

Before any cleaning, the scale must be disconnected from the power supply. The appliance may only be cleaned with a damp cloth. Under no circumstances should water be allowed to enter the unit. Disinfectant may only be used on the display's membrane keypad.

Clean the scale as required with a household cleaner or commercially available disinfectant. Spraying the unit and the connection plug is not permitted.

### 4.6 Maintenance and service

The metrological control was carried out during manufacture by initial calibration. Further regular metrological checks (recalibration) are to be carried out by the responsible calibration authorities in accordance with the respective national rules.

This scale may only be opened and repaired by trained service partners authorised by Soehnle Industrial Solutions. If the scale does not operate as intended, there is a suspicion of damage.

The scale must then be taken to a service partner authorised by Soehnle Industrial Solutions. Only original spare parts may be used for repairs by an authorised service partner. The original parts are described in the service documentation with order number.

### 4.7 Guarantee / Warranty / Liability

Insofar as there is a defect in the delivered item for which the manufacturer is responsible, the manufacturer is entitled to choose either to remedy the defect or to supply a replacement. Replaced parts shall become the property of the manufacturer. If the remedy of the defect of the replacement delivery fails, the statutory provisions shall apply.

The warranty period is **24 months** and starts on the day of the first purchase of the product. **Please keep the invoice as proof.** In case of service, please contact your dealer or the manufacturer's customer service.

**In particular, no warranty is given for damage resulting from the following reasons:**

Unsuitable, improper storage or use, faulty assembly or commissioning by the purchaser or by third parties, natural wear and tear, modification or tampering, faulty or negligent handling, in particular excessive stress, chemical, electrochemical,

## 4. General information

---

electrical influences or moisture, insofar as these are not due to a fault of the manufacturer. If operational, climatic or other influences lead to a significant change in the conditions or the condition of the material, the warranty for the flawless overall function of the units shall lapse. The warranty period for wearing parts (e.g. batteries) is 6 months. **Keep the original packaging for possible return transport!**

### 4.8 Battery and accumulator disposal

Batteries and rechargeable batteries containing harmful substances are marked with the symbol of a crossed-out wheeled bin and must not be disposed of with household waste. As a consumer, you are legally obliged to return used batteries and rechargeable batteries. You can hand in your old batteries and accumulators as hazardous waste at the public collection points in your municipality or wherever batteries of the type in question are sold. You will find these signs on batteries containing harmful substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury. Rechargeable battery packs and lithium batteries may only be disposed of when discharged. If the batteries are not fully discharged, the terminals must be insulated to prevent short circuits, e.g. by covering them with adhesive tape.


### 4.9 Disposing of the scale

According to the current state of knowledge, the device does not contain any special environmentally hazardous substances. This product is not to be treated as normal waste, but must be handed in at a collection point for the recycling of electrical and electronic equipment. For further information, please contact your local authority, the municipal waste disposal services or the company from which you purchased the product.

# 5. Labelling

The illustration shows an example of a type plate of the platform scale 67xx.

The illustration shows an example of a type plate of the platform scale 68xx.

**CE** EC conformity mark

 Symbol protection class II

**REF** Item number of the product

 Manufacturer of the product

**Max** Maximum load of the weighing range

**d=** Digit increment for not approvable scales

**SN** Serial number of the scale

## 6. Assembly and installation - chair scale


---

### 6.1 Unpacking the scale

1. Remove the box from the outer packaging.


2. Remove the polystyrene cover and take out the seat.

3. Put the seat aside for assembly.

4. Remove the front and rear polystyrene cover and the polystyrene cover of the handle.

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**DO NOT DROP THE CHAIR.**

Move the chair to the area designated for assembly.

Put the polystyrene parts back in the box.

# 6. Assembly and installation - chair scale

## 6.2 Mounting the scale

Tools needed for assembly:  
5 mm Allen key

1. loosen both screws for fastening the rear handle frame (marked with a circle) with the Allen key supplied.

**Caution: Do not remove screws.**

2. Remove the polystyrene from the main handle and turn the handle so that the indicator is facing upwards.

3. Insert the main handle into the frame and fix it with the screws.

**Caution: Insert the load cell cable into the grooves provided to avoid damaging the cable.**




4. Turn the armrests until they are parallel to the seat.

**Attention: The armrests can be folded upwards to make it easier to climb onto the scale. Do not lean on the armrests when standing in front of the chair. Unnecessary load can damage the folding mechanism of the armrests. Lock the parking brakes before a patient sits on the chair.**

5. place the seat on the lower part and fix it with four hexagon socket screws M6 x 16.

**Attention: Allen key and the four screws are in the box.**


6. Remove the transport lock (cardboard) from the charging cell.


# 7. Assembly and installation - wheelchair scales

## 7.1 Mounting the scale

Please lift the unpacked scale with the holding/transport handle and screw in the four supplied leveling feet on the underside of the scale.


Using the 4 feet/stand elements, position the scale in a stable, level alignment with the ground.


To set up the scale, operate the levers of the joint locks.  
Seen from the back of the scale:

OPEN the right-hand joint by turning the lever downwards.  
OPEN the left joint by turning the lever upwards.

Fold the railing into the vertical position.

Closing the joint locks:  
CLOSE the right joint by turning the lever upwards.  
CLOSE the left joint by turning the lever downwards.

# 7. Assembly and installation - wheelchair scales

---

## 7.2 Setting up the scale

Place the scale on a flat and hard surface without any unevenness.

### **Attention!**

The weighing platform must touch the floor with all 4 feet, must not wobble and must not touch any other objects.

Set the optimum position of the weighing platform by means of the adjustable foot adjustable foot straps and the integrated spirit level.

### **Hint:**

During weighing, please make sure that the person being weighed does not touch the display and thus affect the result.

## 7.3 Transporting the scale

### **Attention!**

Make sure that no objects are on the scale during transport!

To avoid injury and damage, do not drop or hit the scale during transport.

The scale is equipped with small wheels which are suitable for and achieve transport over shorter distances. Please use this aid.

Remove the power supply before transport.

Fold down the support railing and fix it in the folded position (fix it by means of the turning levers).

Use the extendable handle for transport.

Roll the scale to the desired location. Be careful not to drop the scale during this process as the electronics could be affected. The handle also makes it easier for you to manoeuvre, e.g. in curves.

## 8. Assembly and installation - person scales

---

### 8.1 Setting up the scale

Place the scale on a flat and hard surface without any unevenness.

**Attention!**

The weighing platform must touch the floor with all 4 feet, must not wobble and must not touch any other objects.



Set the optimum position of the weighing platform by means of the adjustable foot adjustable foot straps and the integrated spirit level.

**Hint:**

During weighing, please make sure that the person being weighed does not touch the display and thus affect the result.

# 9. Product description

## 9.1 Display

# 9. Product description

## 9.2 Operating elements

Button	Description	Function
	ON/OFF [Ein/Aus]	Switches the scale on or off.
	PRINT [Drucken] kg/lb	Currently no function
	Zero [Null]	Resets the scale to zero. The reset to zero within three seconds. This function is only possible when the weight is displayed stably. Zeroing up to 2% of the full weight.
	Hold / Release [Halten / Entlasten]	If „HOLD“ is shown in the display, pressing this key displays the current weight value together with the „HOLD“ indicator. If „HOLD“ is not shown in the display and the scale is unloaded, the „HOLD“ display goes out. This key is not active in BMI mode.
	BMI	Call up the BMI mode. Only possible if a weight is shown on the display and this function has been activated in the settings.
	UP / DOWN [Nach oben / unten]	To set the value/number that flashes.
	CLEAR [LÖSCHEN]	When the height is displayed in BMI mode, the value is reset to the default value (170.0 cm). If this button is pressed when the calculated final BMI value is displayed, the scale returns to weighing mode.
	ENTER [EINGABE]	Confirmation of the last parameter entered and call up of the next step. Pressing and holding this key while the scale is being switched on the ID display is called up (pre-parameter mode).
	TARE [TARA]	If the weight is below the maximum value (ASTART value), the manual tare setting is called up. If the weight is above the maximum value, the tare function is called up. During the tare function the indicators NET [Net] & GROSS [Gross] light up.

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## 10.2 Weighing

Stand or sit on the scale. The weight appears in the display.

With the 6702 chair scale, the patient's arms and legs must rest on the arm and leg rests.

## 10.3 Tare

### **Manual tare**

Place the additional weight (weight must be >1kg)) on the weighing platform and then press the tare button. The display flashes and then the display shows „0.0“ and NET (net).

### **Manual tare input**

With the manual tare input you can enter a fixed tare value manually.

With the scale unloaded, press the tare button.

The TARE default value (15 kg) is displayed and the zero flashes.

Set the value with the UP and DOWN arrow keys.

Confirm the entered value with the Enter key. NET (net) is now shown in the display instead of GROSS (gross).

### **Delete tare**

To delete the tare value, press the ZERO key.

## 10.4 Zeros

Press the ZERO button to correct small deviations from the zero point, e.g. due to soiling of the scale.

# 10. Basic functions

---

## 10.5 BMI function

The person steps on the scale while the Zero flag is displayed.  
The weight is displayed.

Press the BMI button. The messages „BMI“ and „CM“ appear on the display and a default value of 170.0 cm (i.e. 170.0) flashes.

The height value is adjusted with the UP/DOWN keys. Each time you press, the value changes by 0.5 cm. If you keep the key pressed, the respective values are increased or decreased continuously.

Pressing ENTER switches to the BMI value.

When ENTER is pressed, the value is accepted. The BMI value appears on the display. Press CLEAR to return to the weighing mode. This switches off the BMI function.

## 10.6 Switching off the scale

Press the ON/OFF button to switch off the scale.

Hint:

Unloaded, the scale switches off automatically after approx. 120 seconds.

# 11. Communication protocols

---

An EDP connection is currently not permissible in accordance with the medical and regulations is currently not permissible.

## 12. Malfunctions - causes and elimination

Read the following troubleshooting instructions before contacting service personnel.

<b>Error</b>	<b>Possible cause</b>	<b>Remedial action</b>
The scale cannot be switched on when using the rechargeable battery. does not switch on when the rechargeable battery is used.	Battery discharged.	Connect the scale to the power source.
The scale cannot be switched on when using the rechargeable battery. does not switch on when the rechargeable battery is used.	<ol style="list-style-type: none"> <li>1. Battery discharged</li> <li>2. Mains socket defective</li> <li>3. Power supply defective</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace batteries</li> <li>2. Use another mains socket</li> <li>3. Replace mains adapter</li> </ol>
Displayed weight questionable or scale does not reset to zero.	<ol style="list-style-type: none"> <li>1. Weight indication is falsified by falsified by the object</li> <li>2. The weighing pan is not correctly positioned</li> <li>3. Before weighing, the display did not show 0.0 „lb/kg was not shown in the display</li> <li>4. Scale is not stable</li> <li>5. Scale is no longer correctly calibrated</li> </ol>	<ol style="list-style-type: none"> <li>1. remove all objects / infant from the weighing pan</li> <li>2. arrange the weighing pan correctly</li> <li>3. remove the infant from the scale, Set the scale to zero and weigh again</li> <li>4. place the scale on a stable surface and repeat weighing procedure</li> <li>5. check weight with object of known with a known weight</li> </ol>
The message „STOP“ appears in the display.	Load exceeds the capacity of the scale.	Remove excessive weight and Observe maximum load values.
Message „---“ is shown in the is shown in the display.	Scale is in underload condition.	Check that the weighing pan on the base.
LO Bat“ is shown in the display. is shown in the display.	Rechargeable battery is weak.	Recharge the battery Recharge.

### Error list

<b>Error</b>	<b>Definition</b>	<b>Check / Replace</b>
SAT	A/D saturation	Check the charging cell - CALL IN A TECHNICIAN
RTC	Reading error of the real-time clock	Battery of real-time clock empty or circuit board defective - CALL IN A TECHNICIAN

In case of non-elimination of the errors or other error messages, please notify your service partner.

## 13. Notes on electromagnetic compatibility

<b>Guidelines and manufacturer's declaration - Electromagnetic emissions</b>		
The 6702 and 6708 series scales are intended for use in the ELECTRICAL ENVIRONMENT specified below. The customer or user of the 6702 and 6708 scales should ensure that it is used in such an environment.		
<b>Emission measurement</b>	<b>Übereinstimmung</b>	<b>Elektromagnetische Umgebung - Leitlinie</b>
HF emissions according to CISPR 11	Group 1	The 6702 and 6708 scales use RF energy only for its internal FUNCTION. Therefore, RF emission is very low and is unlikely to interfere with neighbouring electronic equipment.
HF-Transmissions according to CISPR 11	Class B	Scales 6702 and 6708 are intended for use in all establishments, including residential and those directly connected to a PUBLIC SUPPLY network that also supplies buildings used for residential purposes.
Harmonics according to IEC 61000-3-2	Class A	Scales 6702 and 6708 are intended for use in all establishments, including residential and those directly connected to a PUBLIC SUPPLY network that also supplies buildings used for residential purposes.
Voltage fluctuations / flicker according to IEC 61000-3-3	fulfils	Scales 6702 and 6708 are intended for use in all establishments, including residential and those directly connected to a PUBLIC SUPPLY network that also supplies buildings used for residential purposes

The scale is subject to special precautions regarding EMC and must be installed and commissioned in accordance with the EMC instructions contained in the accompanying PAPERS. Portable and mobile RF communications equipment may affect the scale if the distance is too short.

# 13. Notes on electromagnetic compatibility

<b>Guidelines and manufacturer's declaration - Electromagnetic immunity</b>			
The 6702 and 6708 series scales are intended for use in the ELECTRICAL ENVIRONMENT specified below. The customer or the user of the scales types 6702 and 6708 should ensure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - Guidelines</b>
Static electricity discharge (ESD) according to IEC 61000-4-2	+ 6 kV Contact discharge (indirect)  + 8 kV Air discharge	+ 6 kV Contact discharge  + 8 kV Air discharge	Floors should be made of wood or concrete or have ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30 %.
Fast transient electrical disturbances Bursts according to IEC 61000-4-4	+ 2 kV for mains lines	+ 2 kV for mains lines	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Surges according to IEC 61000-4-5	+ 1 kV Voltage phase conductor - phase conductor	+ 1 kV Voltage phase conductor - phase conductor	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Voltage dips, short-time interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	< 5% UT for ½ period (> 95% slump)  40% UT for 5 period (60 % slump)  70% UT for 25 period (30% dip)  < 5% UT for 5 s (> 95 % dip)	< 5% UT for ½ period (> 95% slump)  40% UT for 5 period (60 % slump)  70% UT for 25 period (30% dip)  < 5% UT for 5 s (> 95 % dip)	The quality of the power supply should be that of a typical commercial or hospital environment. If the user of the 6702 and 6708 scales requires continued FUNCTION even when power interruptions occur, it is recommended that the 6702 and 6708 scales be powered from an uninterruptible power supply.
Magnetic field at the supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should correspond to typical values found in business and hospital environments.
<b>NOTE: UT is the AC mains voltage before the test level is applied.</b>			

# 13. Notes on electromagnetic compatibility

## Guidelines and manufacturer's declaration - Electromagnetic immunity

The 6702 and 6708 series scales are intended for use in the ELECTRICAL ENVIRONMENT specified below. The customer or the user of the scales types 6702 and 6708 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
<p>Conducted RF disturbances according to IEC 61000-4-6</p> <p>Radiated RF disturbances according to IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 26 MHz to 2.5 GHz</p>	<p>Portable and mobile radios shall not be used at a distance from the 6702 and 6708 scales, including lines, less than the recommended separation distance calculated using the equation appropriate for the frequency of transmission.</p> <p>Recommended protective distance:</p> <p><math>d = 0.4 \sqrt{P}</math></p> <p style="padding-left: 150px;">for 80 MHz to 800 MHz</p> <p><math>d = 0.7 \sqrt{P}</math></p> <p style="padding-left: 150px;">for 800 MHz to 2.7 GHz</p> <p>where P is the nominal power of the transmitter in watts (W) as specified by the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>The field strength of stationary radio transmitters is lower than the compliance level at all frequencies according to an on-site survey.</p> <p>Interference is possible in the vicinity of equipment bearing the following pictogram.</p> <div style="text-align: center;">  </div>

NOTE 1: For 26 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not be applicable in all cases. The propagation of electromagnetic magnitudes is affected by absorption and reflection from buildings, objects and people.

- a) The field strength of stationary transmitters, such as base stations of radio telephones and land mobile radios, amateur radio station, AM and FM radio and television transmitters cannot be predicted theoretically with accuracy. To determine the electromagnetic environment with respect to stationary transmitters, a site study should be considered. If the measured field strength at the site where the unit is used exceeds the above compliance levels, the unit should be observed to demonstrate intended operation. If unusual performance characteristics are observed, additional measures may be required, such as modification or relocation of the unit.
- b) Over the frequency range 150 kHz to 80 MHz, the field strength should be less than 3 V/m..

# 13. Notes on electromagnetic compatibility

## Recommended protective distances between portable and mobile HF telecommunications equipment and the scales of type series 6702 and 6708

The 6702 and 6708 series scales are intended for use in the ELECTRONIC ENVIRONMENT specified below. The customer or the user of the scales types 6702 and 6708 can help to avoid electromagnetic interference by maintaining the minimum distance between portable and mobile RF telecommunication devices (transmitters) and the scales types 6702 and 6708 - depending on the output line of the communication device, as specified below.

Nominal power of the transmitter W	Protective distance depending on transmission frequency m		
	150 kHz bis 80 MHz  $d = 1,2 \sqrt{P}$	80 MHz bis 800 MHz  $d = 0,35 \sqrt{P}$	800 MHz bis 2,5 GHz  $d = 0,7 \sqrt{P}$
0,01	0,12	0,04	0,07
0,1	0,38	0,11	0,22
1	1,20	0,35	0,70
10	3,79	1,11	2,21
100	4,0	3,50	7,00

For transmitters whose rated power is not specified in the above table, the distance can be determined using the equation associated with the respective column, where P is the rated power of the transmitter in watts (W) as specified by the transmitter manufacturer.

**NOTE 1:**

An additional factor of 10/3 has been used to calculate the recommended separation distance of transmitters in the frequency range of 80MHz to 2.5 GHz to reduce the likelihood that a mobile/portable communication device inadvertently introduced into the patient area will cause interference.

**NOTE 2:**

These guidelines may not be applicable in all cases. The propagation of electromagnetic quantities is affected by absorption and reflection from buildings, objects and people.

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All rights reserved.
**Soehnle Industrial Solutions GmbH**