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Table of content

1. Scale models	4
2. Scope of delivery.....	4
3. Warnhinweise	5
4. Warnings	6
4.1 Technical data.....	6
4.2 Intended use	7
4.3 Classification.....	7
4.4 Safety instructions.....	7
4.5 Cleaning.....	8
4.6 Maintenance and service	9
4.7 Warranty / Guarantee / Liability.....	9
4.8 Battery and accumulator disposal	10
4.9 Disposal of the scale.....	10
5. Labeling	11
6. Assembly and installation	12
6.1 Unpacking the scale	12
6.2 Components of the scale	13
6.3 Assembly of the scale.....	14
6.4 Inserting the batteries.....	16
6.5 Connecting the power supply	16
7. Product description	17
7.1 Display	17
7.2 Operating elements.....	18
8. Operation	19
8.1 Switching on the scale.....	19
8.2 Weighing.....	19
8.3 Zeroing	19
8.4 Measuring the body height	19
8.5 BMI function	20
8.6 Switching off the scale.....	21
9. Communication protocols	21
10. Malfunctions - causes and elimination	22
11. Note on electromagnetic compatibility	23

1. Scale models

Conformity assessed

Not conformity assessed

These operating instructions apply to the following scale models:

- ▶ 6730.01.002 Personal scale with stand
- ▶ 6731.01.002 Personal scale with stand & height rod
- ▶ 6830.01.002 Personal scale with stand
- ▶ 6831.01.002 Personal scale with stand & height rod

2. Scope of delivery

Personal scale with stand 6730:

Article	Item number
Personal scale with stand	6730.01.002
Power supply	618.020.102
Instruction manual	470.051.220

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Personal scale with stand & height rod 6731:

Article	Item number
Personal scale with stand & height rod	6731.01.002
Power supply	618.020.102
Instruction manual	470.051.220

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3. Warnings

- ▶ In case of malfunction, contact the dealer or the manufacturer. Unauthorized modifications or repairs may damage your scale and void the manufacturer's warranty.
- ▶ The display unit and scale must not get wet. Liquids (e.g. water) can cause damage. Use a dry cloth, e.g. a towel, to dry the display.
- ▶ Disconnect the power supply to this device before any installation, cleaning or maintenance. Otherwise, the device could be damaged.
- ▶ If the scale is not to be used for a long period of time, please remove the power supply.
- ▶ Avoid stacking materials on the display or loading the display with weights. This can lead to damage.
- ▶ Place the scale on a firm, stable and level surface so that accurate measurement results can be ensured. If the surface is soft or inclined, the measurement results will not be representative.
- ▶ Do not connect the display to unstable power sources.
- ▶ Use only the original equipment. The use of other makes may cause damage to the scale.
- ▶ Simultaneous contact between the interface and the patient is not permitted.

4. General notes

4.1 Technical data

	Personal scale with stand 6730:	Personal scale with stand 6830:
Scale type	Single range scale	Single range scale
Maximum load	250 kg	250 kg
Minimum load	2 kg	2 kg
Digit increment	100 g	100 g
Dimensions (W x D x H)	375 x 375 x 1.330 mm	375 x 375 x 1.330 mm
Product weight	10,5 kg	10,5 kg
Power supply (6730)	Adapter 230V~/120V~ 60/50Hz-9VDC 500mA Use only power supplies that are approved for medical applications applications and comply with the EN 60601 standard.	Adapter 230V~/120V~ 60/50Hz-9VDC 500mA
Calibration class	Calibration class III, MPG class 1 m	
Working temperature	0 °C bis + 40 °C	0 °C bis + 40 °C
Storage temperature	0 °C bis + 50 °C	0 °C bis + 50 °C
Humidity	85 % (non-condensing)	85 % (non-condensing)
Air pressure	Use at atmospheric pressure	Use at atmospheric pressure

	Personal scale with stand & height rod 6731:	Personal scale with stand & height rod 6831:
Scale type	Single range scale	Single range scale
Maximum load	250 kg	250 kg
Minimum load	2 kg	2 kg
Digit step	100 g	100 g
Measuring range	75 -205 cm	75 -205 cm
Dimensions (W x D x H)	375 x 375 x 1.330 mm	375 x 375 x 1.330 mm
Product weight	10,5 kg	10,5 kg
Power supply	Adapter 230V~/120V~ 60/50Hz-9VDC 500mA Use only power supplies that are approved for medical applications applications and comply with the EN 60601 standard.	Adapter 230V~/120V~ 60/50Hz-9VDC 500mA
Calibration class	Calibration class III, MPG class 1 m	
Working temperature	0 °C bis + 40 °C	0 °C bis + 40 °C
Storage temperature	0 °C bis + 50 °C	0 °C bis + 50 °C

4. General notes

Humidity	85 % (non-condensing)	85 % (non-condensing)
Air pressure	Use at atmospheric pressure	Use at atmospheric pressure

4.2 Intended use

These scales are used to determine the weight of persons and may only be used for this purpose. Any other use is prohibited. These scales are intended for use in legal metrology. The scales comply with the provisions of the applicable European directives and regulations.

The manufacturer's medical product advisor must be notified in writing of any faults which could endanger the patient and of any faults which could lead to the falsification of measurement results. These instructions for use are an integral part of the equipment. Precise observance of these instructions is a prerequisite for the intended use and correct operation of the devices.

The devices may only be operated by qualified personnel. Operation by the patient is not intended.

Please ensure that the scales are not subjected to permanent loads. In addition, shocks and overloading of the scales in excess of the specified maximum load must be avoided at all costs. This may damage the scales.

4.3 Classification

Electrical protection class II (protective insulation, no protective earth connection).

Application part (6730 & 6731)

The scale is a Class 1 medical device with a measuring function. The medical device is a type B applied part and is intended for direct contact by the patient. The leakage currents correspond to the classification of type B applied parts

4. General notes

4.4 Safety instructions

These scales are designed for medical use. The user must be familiar with the operation of the scale. Before operating the scale, please read the information in the operating instructions carefully. It contains important instructions for installation, intended use and maintenance of the device. The manufacturer is not liable if the following instructions are not observed. Improper installation will void the warranty.

- ▶ The scale is not explosion-proof and must not be used in critical or potentially explosive atmospheres.
- ▶ When using electrical components under increased safety requirements, the corresponding regulations must be observed.
- ▶ Electrical connection conditions must match the values printed on the power supply unit.
- ▶ The scale is designed for operation in buildings. Observe the permissible ambient conditions for operation (see Technical data). The scale meets the requirements for electromagnetic compatibility. Exceeding the maximum values specified in the standards must be avoided.
- ▶ The scale must be set up in such a way that disconnection from the mains can be easily accomplished at any time. The accessibility of the socket should be ensured at all times.
- ▶ When moving the scale, lift the scale. The scale must never be pushed. Then check the leveling and adjust it to the new situation if necessary.
- ▶ The scale is a measuring instrument. Drafts, vibrations, rapid temperature changes and solar radiation can affect the weighing result. High humidity, vapors, aggressive liquids and heavy soiling must be avoided.
- ▶ This device is interference suppressed in accordance with the applicable EC directive. However, under extreme electrostatic and electromagnetic influences, e.g. when operating a radio or cell phone in the immediate vicinity of the device, the display value may be affected. After the interference has ceased, the product can be used again as intended; if necessary, it must be switched on again.

4.5 Cleaning

Before each cleaning, the scale must be disconnected from the power supply. The

4. General notes

device may only be cleaned with a damp cloth. Under no circumstances may water penetrate the device. Disinfectant may only be used on the membrane keypad of the display.

Clean the scale as required with a household cleaner or commercially available disinfectant. Spraying the device and the connector plug is not permitted.

The metrological control was carried out during the production by the initial verification carried out. Further regular metrological checks (recalibration) must be carried out by the responsible verification authorities in accordance with the respective national regulations.

4.6 Maintenance and service

This scale may only be opened and repaired by trained service partners authorized by Soehnle Industrial Solutions. If the scale does not operate as intended, there is a suspicion of damage.

The scale must then be taken to a service partner authorized by Soehnle Industrial Solutions. Only original spare parts may be used for repairs by an authorized service partner. The original parts are described in the service documentation with order number.

4.7 Warranty / Guarantee / Liability

Insofar as there is a defect in the delivered item for which the manufacturer is responsible, the manufacturer is entitled to choose either to remedy the defect or to deliver a replacement. Replaced parts become the property of the manufacturer. If the remedy of the defect of the replacement delivery fails, the statutory provisions shall apply.

The warranty period is 24 months and starts on the day of the first purchase of the product. Please keep the invoice as proof. In case of service, please contact your dealer or the manufacturer's customer service.

In particular, no warranty is given for damage resulting from the following causes: Unsuitable, improper storage or use, faulty assembly or commissioning by the purchaser or by third parties, natural wear and tear, modification or tampering, faulty or negligent handling, in particular excessive stress, chemical, electrochemical, electrical influences or moisture, unless these are the fault of the manufacturer. If operational, climatic or other influences lead to a significant change in the conditions or the condition of the material, the warranty for the proper overall function of the equipment shall become void. The warranty period for wearing parts

4. General notes

(e.g. rechargeable batteries) is 6 months.

Keep the original packaging for possible return transport!

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
Rechargeable battery packs and lithium batteries may only be disposed of when discharged. If the batteries are not fully discharged, the terminals must be insulated to prevent short circuits, e.g. by covering them with adhesive tape.

4.9 Disposing of the scale

According to the current state of knowledge, the device does not contain any special environmentally hazardous substances. This product is not to be treated as normal waste, but must be taken to an acceptance point for the recycling of electrical and electronic equipment. Further information can be obtained from your local authority, the municipal waste disposal companies or the company from which you purchased the product.

5. Marking

The illustration shows an example of a type plate of the platform scale 67xx.

The illustration shows an example of a type plate of the platform scale 68xx.

CE EC conformity mark

 Symbol protection class II

REF Item number of the product

 Manufacturer of the product

Max Maximum load of the weighing range

d= Digit increment for not approvable scales


SN Serial number of the scale

6. Assembly and mounting

6.1 Unpacking the scale

Visually inspect the scale and verify that all components are included in the carton.

The carton should contain the scale, these instructions, and a box labeled „PARTS INSIDE“.

The box contains the following parts:

- Two (2) wheels.
- One (1) wheel axle.
- Five (5) M5 x 35 flat head Phillips screws.
- Two (2) M4 x 16 flat head Phillips screws.

- Power supply 230/120 VAC, 50/60 Hz - 9 VDC, 0.5 A
- 6 AA batteries

6. Assembly and mounting

6.2 Components of the scale

Remove the scale components from the box and place them clearly arranged on a free work surface.

Weighing platform

Upper stand with display

Length gauge


Lower stand

Tools needed for assembly:
Phillips screwdriver, size 2

6. Assembly and mounting

6.3 Mounting the scale


1. Insert the lower stand into the weighing platform, using the bolts as guides. the bolts as guides.


Be careful that the load cell cable does NOT get between these parts.

Make sure that there are two screw holes on the bottom stand that should be facing up when reassembling.

2. Insert the four M5 screws and tighten them with a Phillips screwdriver.

3. Remove the two wheels and the axle from the box labeled „PARTS INSIDE“. Position the wheels on the axle, mount them in the groove, and lock them in place with the M5 screw using a Phillips screwdriver.


6. Assembly and mounting

4. Insert the upper stand into the lower stand.


Be careful not to damage the load cell cable.

Insert the edge of the length gauge into the appropriate rail.

Attach the upper instrument assembly by screwing the two M4 screws into the back of the instrument back of the unit.

5. Connect the upper length gauge to the lower length gauge.

6. Assembly and mounting

6.4 Inserting the batteries

1. Open the battery compartment cover.

2. insert the supplied batteries.


3. Close the battery compartment cover.



6.5 Connecting the power supply unit

The power supply is connected to the right socket marked with the specification for the input voltage. The power supply unit is then plugged into the socket.

7. Product description

7.1 Display

7. Product description

7.2 Operating elements

Button	Designation	Function
	ON/OFF [Ein/Aus]	Turns the scale on or off.
	PRINT [Drucken] kg/lb	Currently no function.
	Zero [Null]	Resets the scale to zero. The reset to zero within three seconds. This function is only possible if the weight is displayed stable. Zeroing up to 2% of the full weight.
	Hold / Release [Halten / Entlasten]	<p>If „HOLD“ is shown in the display, the current weight value appears together with the „HOLD“ display by pressing this key.</p> <p>If „HOLD“ is not shown in the display and the scale is unloaded, the „HOLD“ display goes out.</p> <p>This key is not active in the BMI mode.</p>
	BMI	Calling up the BMI mode. Only possible if a weight is shown on the display and this function has been activated in the settings.
	UP / DOWN [Nach oben / unten]	To set the value/number that flashes.
	CLEAR [LÖSCHEN]	<p>When the height is displayed in BMI mode, the value is reset to the default value (170.0 cm).</p> <p>If this key is pressed when the calculated final BMI value is displayed, the scale returns to the weighing mode.</p>
	ENTER [EINGABE]	<p>Confirmation of the last parameter entered and call of the next step.</p> <p>By keeping this key pressed during the power-on procedure of the scale the ID display is called up (pre-parameter mode).</p>

8. Basic functions

8.1 Switching on the scale

With the scale unloaded, press the ON/OFF key. All symbols and digits first appear on the display. Then „00.0“ appears in the display.

The scale is now ready for weighing.

8.2 Weighing

Ask the patient to step on the scale carefully.
The weight appears automatically after the scale is loaded.

When weighing, both legs of the patient must be on the weighing platform of the scale.


8.3 Zeros

Press the ZERO key to correct small deviations from the zero point, e.g. due to contamination of the scale.

8.4 Measurement of body size

People between 1.09 m and 2.00 m tall:

Unfold the extension arm of the length gauge and carefully push it upwards. Position the patient in the center of the scale. Then slowly push the unfolded length measuring rod down again until the extension arm rests on the head of the person to be measured. The attached scale can now be used to read off the body height in cm.


8. Basic functions

8.4 Measurement of body size

People between 70 cm and 1.09 m tall:

Unfold the extension arm of the length gauge. Position the patient in the center of the scale. Then slowly push the unfolded length measuring rod down until the extension arm rests on the head of the person to be measured. The attached scale can now be used to read off the body height in cm.

8.5 BMI function

The person steps on the scale while the Zero flag is displayed.
The weight is displayed.

Press the BMI key. The messages „BMI“ and „CM“ appear on the display and a default value of 170.0 cm (i.e. 170.0) flashes.

The height value is adjusted with the UP/DOWN (ARROW UP/ARROW DOWN) keys. Each time the key is pressed, the value changes by 0.5 cm. If you keep the key pressed, the respective values are continuously increased or decreased.

Pressing ENTER switches to the BMI value.

When ENTER is pressed, the value is accepted. The BMI value appears on the display. Press CLEAR to return to the weighing mode. This will turn off the BMI function.

8. Basic functions

8.6 Switching off the scale

Press the ON/OFF key to turn the scale off.

Notice:

Unloaded, the scale switches off automatically after approx. 120 seconds.

9. Communication protocols

An EDP connection is currently not permissible in accordance with the medical and regulations is currently not permitted.

10. Malfunctions - causes and elimination

Read the following troubleshooting instructions before contacting service personnel.

Error	Possible cause	Remedial action
The scale cannot be switched on when using the rechargeable battery. Battery cannot be switched on.	Battery discharged.	Connect the scale to the power source.
The scale cannot be switched on when using the rechargeable battery. Does not turn on when the rechargeable battery is used.	<ol style="list-style-type: none"> 1. Battery discharged 2. Mains socket defective 3. Power supply defective 	<ol style="list-style-type: none"> 1. Replace batteries 2. Use other power outlet 3. Replace power supply
Displayed weight questionable or scale does not reset to zero.	<ol style="list-style-type: none"> 1. Weight display is falsified by object 2. The weighing pan is not correctly positioned 3. Before weighing, the display did not show 0.0 „lb/kg was not displayed before weighing 4. Scale is not stable 5. Scale no longer correctly calibrated 	<ol style="list-style-type: none"> 1. All objects / infant from the weighing pan 2. Arrange weighing pan correctly 3. Remove infant from weighing pan, Set scale to zero and weigh again 4. Place scale on stable surface and repeat weighing procedure 5. Check weight with object of with known weight
„STOP“ message is shown in the display.	Load exceeds the capacity of the scale.	Remove excessive weight and Observe maximum load values.
Message „---“ is shown in the display.	Scale is in underload condition	Check that the weighing pan on the base.
"LO Bat" is shown in the display.	Rechargeable battery is weak.	Recharge the battery recharge.

Error list

Error	Definition	Check / Replace
SAT	A/D Saturation	Check charging cell - CALL IN TECHNICIAN
RTC	Real time clock reading error	Battery of the real-time clock empty or circuit board defective - CALL IN A TECHNICIAN

If the errors cannot be eliminated or if other error messages are displayed, please contact your your service partner.

11. Notes on electromagnetic compatibility

Guidelines and manufacturer's declaration - Electromagnetic emissions		
The scales of the 6730 and 6731 series are intended for operation in the ELECTRICAL ENVIRONMENT specified below. The customer or the user of the scales types 6730, 6731 should ensure that it is used in such an environment.		
Emittance measurement	Match	Electromagnetic environment - Guideline
RF emissions according to CISPR 11	Group 1	The 6730 and 6731 scales use RF energy only for its internal FUNCTION. Therefore, RF emission is very low and is unlikely to interfere with nearby electronic equipment.
RF emissions according to CISPR 11	Class B	The 6730 and 6731 scales are intended for use in all facilities, including residential and those directly connected to a PUBLIC POWER SUPPLY that also supplies buildings used for residential purposes.
Harmonics according to IEC 61000-3-2	Class A	Scales 6730 and 6731 are intended for use in all facilities, including residential and those directly connected to a PUBLIC SUPPLY NETWORK that also serves buildings used for residential purposes.
Voltage fluctuations / flicker according to IEC 61000-3-3	Fulfills	The 6730 and 6731 scales are intended for use in all establishments, including residential, and those directly connected to a PUBLIC SUPPLY NETWORK that also supplies buildings used for residential purposes.

The scale is subject to special precautions regarding EMC and must be installed and put into operation in accordance with the EMC instructions contained in the accompanying PAPERS. Portable and mobile RF communications equipment can affect the scale if the distance is too short.

11. Notes on electromagnetic compatibility

Guidelines and manufacturer's declaration - Electromagnetic immunity			
The scales of type series 6730 and 6731 are intended for operation in the ELECTRICAL ENVIRONMENT specified below. The customer or the user of the scales types 6730 and 6731 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	+ 6 kV contact discharge (indirect) + 8 kV Air discharge	+ 6 kV contact discharge + 8 kV Air discharge	Floors should be made of wood or concrete, or have ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Fast transient electrical disturbance magnitudes/bursts according to IEC 61000-4-4	+ 2 kV for power lines	+ 2 kV for power lines	The quality of the supply voltage should be that of a typical business or hospital environment.
Surges according to IEC 61000-4-5	+ 1 kV Voltage outer conductor - outer conductor	+ 1 kV Voltage outer conductor - outer conductor	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short-time interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	< 5% UT for ½ period (> 95% slump) 40% UT for 5 period (60% slump) 70% UT for 25 period (30% dip) < 5% UT for 5 s (> 95% dip)	< 5% UT for ½ period (> 95% slump) 40% UT for 5 period (60% slump) 70% UT for 25 period (30% dip) < 5% UT for 5 s (> 95% dip)	The quality of the supply voltage should be that of a typical commercial or hospital environment. If the user of the 6730 and 6731 scales requires continued FUNCTION even when power interruptions occur, it is recommended that the 6730 and 6731 scales be powered from an uninterruptible power supply.
Magnetic field at supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should correspond to typical values found in business and hospital environments.
NOTE: UT is the AC line voltage before the test level is applied.			

11. Notes on electromagnetic compatibility

Guidelines and manufacturer's declaration - Electromagnetic immunity

The scales of type series 6730 and 6731 are intended for operation in the ELECTRICAL ENVIRONMENT specified below. The customer or the user of the scales types 6730 and 6731 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
<p>Conducted RF disturbances according to IEC 61000-4-6</p> <p>Radiated RF disturbances according to IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 26 MHz to 2.5 GHz</p>	<p>Portable and mobile radios will not be used at a distance from 6730 and 6731 scales, including lines, less than the recommended protective distance calculated using the equation appropriate for the transmit frequency.</p> <p>Recommended guard distance:</p> <p>$d = 0.4 \sqrt{P}$</p> <p style="padding-left: 100px;">for 80 MHz to 800 MHz</p> <p>$d = 0.7 \sqrt{P}$</p> <p style="padding-left: 100px;">for 800 MHz to 2.7 GHz</p> <p>where P is the nominal power of the transmitter in watts (W) as specified by the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>The field strength of stationary radio transmitters is less than the compliance level at all frequencies according to an on-site survey.</p> <p>Interference is possible in the vicinity of equipment bearing the following pictogram.</p> <div style="text-align: center;">  </div>

NOTE 1: At 26 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not be applicable in all cases. The propagation of electromagnetic is affected by absorption and reflection from buildings, objects and people.

a) The field strength of stationary transmitters, such as base stations of radio telephones and land mobile radios, amateur radio station, AM and FM radio and television transmitters cannot be predicted theoretically with accuracy. To determine the electromagnetic environment with respect to stationary transmitters, a site study should be considered. If the measured field strength at the site where the equipment is used exceeds the above compliance levels, the equipment should be observed to demonstrate intended operation. If unusual performance characteristics are observed, additional measures may be required, such as modification or relocation of the device.

b) Over the frequency range 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

11. Notes on electromagnetic compatibility

Recommended protective distances between portable and mobile RF Telecommunications Equipment and the 6730 and 6731 Series Scales

The scales of type series 6730 and 6731 are intended for operation in the ELECTRICAL ENVIRONMENT specified below. The customer or the user of the scales types 6730 and 6731 can help to avoid electromagnetic interference by maintaining the minimum distance between portable and mobile RF telecommunication devices (transmitters) and the scales types 6730 and 6731 - depending on the output line of the communication device, as indicated below.

Nominal power of the transmitter W	Protective distance depending on transmission frequency m		
	150 kHz bis 80 MHz $d = 1,2 \sqrt{P}$	80 MHz bis 800 MHz $d = 0,35 \sqrt{P}$	800 MHz bis 2,5 GHz $d = 0,7 \sqrt{P}$
0,01	0,12	0,04	0,07
0,1	0,38	0,11	0,22
1	1,20	0,35	0,70
10	3,79	1,11	2,21
100	4,0	3,50	7,00

For transmitters whose rated power is not specified in the above table, the distance can be determined using the equation associated with the respective column, where P is the rated power of the transmitter in watts (W) as specified by the transmitter manufacturer.

NOTE 1:

An additional factor of 10/3 was used to calculate the recommended separation distance of transmitters in the frequency range of 80MHz to 2.5 GHz to reduce the likelihood that a mobile/portable communications device inadvertently introduced into the patient area would cause interference.

NOTE 2:

These guidelines may not be applicable in all cases. The propagation of electromagnetic quantities is affected by absorption and reflection from buildings, objects, and people.

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All rights reserved.
Soehnle Industrial Solutions GmbH