

9300N

Electrical Patient Lift & Transfer Chair

The innovative design makes transferring the patient
from the bed to the chair a breeze!



advantages

★ ★ ★ 6 ★ ★ ★



Electrical Lifting
Dual alu column lifting
motors



IP25 waterproof



Easy to open and
close



High load Capacity



Triple security
protection



Medical Silent
Wheels

Products & Details



1

180° Degree easy to open and close



2

IP25 Waterproof

The chair has an IP25 waterproof rating and can be used as a shower chair directly into the bathroom without taking out the battery.



4

Electric up & down

150kg Lfiting Capacity, Lifting stroke 200mm.



5

Foldable handle

The handle can be folded by a quick releaseclamp, reserving more space for daily activities.



6

Triple Safety

The chair has safety belt and Safety buckle and safety latch.



242.5 lbs

The maximum load is 242.5 lbs

3

- Weight 36kg
- Loading Capacity 110kg/242.5lbs
- Actuator 24V/100W×2
- Battery 25.9V/2600mAh
- Waterproof level IP25

