

PUMA AIR ROLLATOR

The **Puma Air** rollator is a high quality outdoor rollator. This rollator has a matte black or red frame and extra large 12" pneumatic tyres (295 mm) for more stability on difficult surfaces. The push handles can easily be adjusted in height without the use of any tools. In addition, the rollator can be folded very easily with a flick of the wrist by pulling the strap on the seat upwards. As standard the Puma Air comes with: a back belt, a cup holder, a crutch holder and a modern, closed shopping bag. With its large pneumatic tyres, the Puma Air is especially suitable for outdoor use, but also for indoor use as it is only 65 cm wide. In total, the outdoor rollator only weighs 9,78 kg and can therefore easily be transported. The rollator is tested and approved for a user weight of up to 150 kg.



Puma Air Rollator	312501+02
Height of the push handles (adjustable)	82 - 99 cm
Recommended user heights	160 - 200 cm
Seat height	59 cm
Total width opened / folded	65 cm / 26 cm
Inner width between push handles	40 cm
Length (depth)	76,4 cm
Weight	9,78 kg
Max. load (rollator / bag)	150 kg / 5 kg
Brakes	2 lockable with parking function
Wheel sizes	Front and back: 295 mm (12") x 45 mm



312560

Accessory:
Tray that is easily attached to the rollator (instead of the back belt).




ACCESSORIES

A small selection of accessories available for our Mobilex rollators. The entire selection can be found in our catalog or on www.mobilex.dk.




Umbrella



One-hand brake



Oxygen holder



Open net bag



Bell



LED lights



High back belt



Tray



Fixed padded seat



Bag with zip



Speed reduction brake

* Only available for the Buffalo, Gepard, Leopard, Lion and Tiger rollator