



Personal scale Conformity assessed

6710

- ▶ Model 6710: Conformity assessed according to Class III, MPG Class 1 m
- ▶ Large tread surface: 375 x 375 mm
- ▶ Easy to clean surface
- ▶ Illuminated LC display, digit height: 20 mm
- ▶ Integrated RS232 and USB interfaces for connecting a PC



Features

- ▶ Model 6710: Conformity assessed according to Class III, MPG Class 1 m
- ▶ Large tread surface: 375 x 375 mm
- ▶ Easy to clean surface
- ▶ Terminal with plastic housing
- ▶ Cable length platform to terminal: 1.8 m
- ▶ Illuminated LC display, digit height: 20 mm
- ▶ Display of weight and BMI
- ▶ Integrated RS232 and USB interfaces for connecting a PC
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 6 x AA batteries and mains adapter operation 100 -240 V, 50 - 60 Hz possible (both included in the scope of delivery)
- ▶ Hold function for convenient reading of the weight value after weighing

Not for delivery to France!

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
6710.01.002	No	100	250
Personal scale			
