



Personal scale with stand conformity assessed

6730

- ▶ Conformity assessed according to class III, MPG class 1 m
- ▶ Backlit LC display, digit height: 20 mm
- ▶ Indication of weight and BMI
- ▶ Battery operation and mains operation
- ▶ Easy transport to the patient



Features

- ▶ Conformity assessed according to class III, MPG class 1 m
- ▶ Large tread: 374 x 374 mm
- ▶ Easy to clean surface
- ▶ Backlit LC display, digit height: 20 mm
- ▶ Automatic calculation of the BMI by entering the body height and the determined weight value
- ▶ Indication of weight and BMI
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 6 x AA batteries and mains operation 100 - 240 V 50 - 60 Hz possible (both included in delivery)
- ▶ Easy transport to the patient due to two castors and a handle on the back of the stand
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Integrated RS232 and USB interfaces for connecting a PC

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
6730.01.002	Yes	100	250
Personal scale with stand			