



Personal scale with stand & height rod conformity assessed

6731

- ▶ Conformity assessed according to class III, MPG class 1 m
- ▶ Integrated height rod
- ▶ Indication of weight, height and BMI
- ▶ Battery operation and mains operation
- ▶ Easy transport to the patient



Features

- ▶ Conformity assessed according to class III, MPG class 1 m
- ▶ Integrated height rod, Unit: 1 mm
- ▶ Large tread: 374 x 374 mm
- ▶ Easy to clean surface
- ▶ Backlit LC display, digit height: 20 mm
- ▶ Automatic calculation of the BMI by entering the body height and the determined weight value
- ▶ Indication of weight, height and BMI
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 6 x AA batteries and mains operation 100 - 240 V, 50 - 60 Hz possible (both included in delivery)
- ▶ Easy transport to the patient due to two castors and a handle on the back of the stand
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Integrated RS232 and USB interfaces for connecting a PC

Products

Order number Variation	Approvable	Division g	Maximum load kg
6731.01.002 Personal scale with stand & height rod	Yes	100	250
