



Baby scale "Neo" conformity assessed

6752

- ▶ Model 6752: Conformity assessed according to class III, MPG class 1 m
- ▶ Very fine division 5 g
- ▶ Easy-to-read LC display with 20 mm high digits and backlight
- ▶ Integrated length scale in the weighing tray
- ▶ Battery and mains operation



Features

- ▶ Model 6752: Conformity assessed according to Class III, MPG Class 1 m
- ▶ Particularly small division 5 g
- ▶ Easy-to-read LC display with 20 mm high digits and backlight
- ▶ Tare function: documents placed on the scale can be tared
- ▶ Tare hand entry possible via the tare button
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Integrated RS232 and USB interfaces for connecting a PC
- ▶ Integrated length scale in the weighing pan (units in 2.5 mm increments)
- ▶ Weighing pan material: Acrylic
- ▶ Units in kg
- ▶ Energy-saving auto-off function
- ▶ Battery and mains operation 100 - 240 V, 50-60 Hz (included in scope of delivery)
- ▶ Battery status indicator

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
6752.01.002	Yes	5	15
Baby scale "Neo" conformity assessed			

Accessories

Order number

Variation

Description

5005.11.001



Baby height rod

2550.03.022

USB data cable, cable length 3 m
