



Handrail scale

6803

- ▶ Stainless steel handrails with padded handles
- ▶ Four leveling feet for levelling and stabilizing the scale
- ▶ Functions: Weighing, taring, hold, BMI
- ▶ Energy-saving auto-off function

Features

- ▶ Powder-coated base frame, platform made of stainless steel
- ▶ Terminal with aluminium housing
- ▶ Large LC display, digit height 25,4 mm
- ▶ Stainless steel handrails with padded handles
- ▶ Four leveling feet for levelling and stabilizing the scale
- ▶ Easy transport to the patient due to two castors
- ▶ Functions: Weighing, taring, hold, BMI
- ▶ Battery level indicator (battery symbol in display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg, lb or lb:oz
- ▶ Serial RS 232 and USB interface
- ▶ Battery operation with four Mignon batteries 1.5 V AA (included in scope of delivery)
- ▶ Power supply operation optional (see accessories)

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
6803.01.001	No	100	300
Handrail scale			

Accessories

Order number

Variation

Description

5908.09.001

Power supply 9V / 2A

2550.03.014

Data cable RS 232 PC/EDV 3 m
