



Wheelchair weigher

6808

- ▶ Platform with integrated drive-on ramps
- ▶ Including foldable railing, white painted
- ▶ Transportable due to two integrated castors and handle
- ▶ Battery and mains operation



Features

- ▶ Platform with integrated drive-on ramps
- ▶ Including foldable railing, white painted
- ▶ Passage width: 760 mm
- ▶ Terminal with plastic housing
- ▶ Backlit LC display, digit height: 20 mm
- ▶ Transportable due to two integrated castors and handle
- ▶ Functions: weighing, taring, hold, BMI
- ▶ Battery level indicator (battery symbol in the display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg
- ▶ Battery and mains operation 100 - 240 V, 50-60 Hz (included in scope of delivery)
- ▶ Integrated RS232 and USB interfaces for connecting a printer or PC

Not for delivery to France!

Products

Order number Variation	Approvable	Division g	Maximum load kg
6808.01.502 Wheelchair weigher	No	100	300
