



Chair weigher

7802

- ▶ Robust construction
- ▶ Simple weighing of persons with a mobility handicap
- ▶ Ergonomically shaped and easy to clean seat
- ▶ Hold and tare function



Features

- ▶ Ergonomically shaped and easy to clean seat
- ▶ Dimensions seat: 41 x 40 cm
- ▶ Seat height: approx. 51 cm
- ▶ Hinged arm- and footrests
- ▶ Mains adapter operation (included in the scope of delivery)
- ▶ Battery operation (6x AA batteries, not included) or optional rechargeable battery operation with battery pack and charging circuit
- ▶ Energiesparende Auto-Off-Funktion
- ▶ Taring function: Placed underlays can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Version .001: two rear swivel castors with rear fixing brakes
- ▶ Version .401: four swivel castors with rear fixing brakes
- ▶ Switchable to kg/lb units

Products

Order number Variation	Approvable	Division g	Maximum load kg	Model
7802.65.001 Chair weigher with Terminal 3810	No	100	200	2 rear swivel castors with rear fixing brakes
7802.65.401 Chair weigher with Terminal 3810	No	100	200	4 swivel castors with rear fixing brakes

Accessories

Order number

Variation

Description

2557.40.05

Battery incl. charging circuit
