



# Handrail scale

7808

- ▶ Low installation height of only 55 mm
- ▶ Including adiposity railing set, painted white
- ▶ Energy-saving auto-off function for battery operation
- ▶ Tare function: placed documents can be tared



## Features

---

- ▶ Large, flat weighing platform made of aluminium, painted white
- ▶ Low installation height of only 55 mm
- ▶ Including adiposity railing set, painted white
- ▶ Includes two castors and carrying handle
- ▶ Optional wall bracket for the terminal
- ▶ Mains adapter operation (included in scope of delivery)
- ▶ Battery operation (6x AA batteries, not included) or optional rechargeable battery operation with battery pack and charging circuit
- ▶ Energy-saving auto-off function for battery operation
- ▶ Tare function: placed documents can be tared
- ▶ Hold function for convenient reading of the weight value after weighing
- ▶ Units in kg/lb

## Products

---

Order number	Approvable	Division	Maximum load
Variation		g	kg
<b>7808.65.701</b>	No	100	300
Handrail scale			

---

## Accessories

---

**Order number**

Variation

**Description**

---

**2503.30.005**

Wall bracket for terminal 3810

---

**2557.40.05**

Battery incl. charging circuit

---