



Personal scale with stand

7835

- ▶ Display of weight and BMI
- ▶ Backlit LC display, digit height 20 mm
- ▶ Mother-child weighing
- ▶ Energy-saving auto-off function
- ▶ Integrated castors for easy transport of the scales



Features

The stand scale is the ideal scale for fitness studios, pharmacies, care or other social institutions. Due to the raised display, the weight value can be read off comfortably at eye level and through illuminated digits. The tread surface is robust and easy to clean. The scale is a special highlight due to its various functionalities. For example, with the help of the hold function, the weight value can also be read off only after weighing. The mother-child weighing function enables the weight of small children who cannot yet stand on the scales themselves to be determined. The scale also calculates the BMI index. After entering your height, you will receive your BMI index immediately during weighing. With the integrated interface the values can be transmitted directly to a PC or printer.

The scale can be operated by a power supply unit or be used with batteries for mobile applications. The automatic switch-off function switches the scale off when not in use, thus ensuring that no unnecessary energy is consumed. Integrated castors also make it easy to transport the scale.

- ▶ Automatic calculation of the BMI by entering the height and the determined weight value
- ▶ Display of weight and BMI
- ▶ Large tread surface: 350 x 350 mm
- ▶ Sturdy construction with aluminium stand
- ▶ Backlit LC display, digit height 20 mm
- ▶ Easy to clean surface
- ▶ Tare function
- ▶ Hold function for convenient reading of the weight value after weighing
- ▶ Mother-child weighing
- ▶ RS232 interface optional, for connecting a printer or PC
- ▶ Unit changeover kg/lb
- ▶ Energy-saving auto-off function
- ▶ Battery operation with six Mignon 1.5 V AA batteries (included in the scope of delivery)
- ▶ Optional mains adapter operation (see accessories)
- ▶ Integrated castors for easy transport of the scales

Products

| Order number Variation | Approvable | Division g | Maximum load kg |
|-----------------------------------|------------|---------------|--------------------|
| 7835.01.001 Stand scale | No | 100 | 200 |
| 7835.01.002 Stand scale | No | 100 | 300 |

Accessories

| Order number Variation | Description |
|---------------------------|--|
| 2557.05.002 | Power supply unit 100 - 240 V, 50-60 Hz |
| 2557.25.003 | UK exchange adapter |
| 2550.03.002 | Data cable (cable length 2 m) |
| 2563.41.005 | Interface RS 232, installed (must be included in the purchase order, no retrofit possible) |