



Examination table with integrated scale

7857

- ▶ Quick and easy 1-button operation of the scale
- ▶ Easy to read LC display, digit height 20 mm
- ▶ Electrical height adjustable table from 40-95 cm with pedal



Features

This veterinary examination table with integrated scale is designed for daily use in veterinary surgeries. As a treatment or operating table it offers the possibility to weigh the animals at the same time. The animal in need of treatment therefore does not have to be relocated to another scale.

- ▶ Quick and easy 1-button operation of the scale
- ▶ Easy to read LC display, digit height 20 mm
- ▶ Hold function: For easy reading of the weight value after weighing
- ▶ Tare function: Covers/paddings can be tared
- ▶ Display made of high-quality and durable plastic
- ▶ Large stainless steel weighing platform (type 304L/18/10)
- ▶ Electrical height adjustable table from 40-95 cm with pedal
- ▶ Construction of a scissors lift table
- ▶ Technical data of the motor: 220 V, 400 Watt
- ▶ Units in kg or lb

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
7857.70.001	No	50	150
Examination table with integrated scale			