


FASTER +ROOF

available colours:



If you have ever felt too hot while enjoying the outdoors, then **FASTER** mobility scooter with roof is just what you need. With a perfect canopy, it will shield you from the heat, and it will make every ride much more pleasant. This model has been especially welcome by all those who love playing golf. Four wheel stability provide great support and the weighted roof keeps the vehicle to the ground while not impacting its other abilities.

The scooter comes **assembled and ready to go!** 

SPECIFICATIONS

OVERVIEW

Max user weight: 160 kg / 25.2 st
Scooter weight: 133 kg / 20.9 st (w/ Lead-Acid Gel battery) / 117 kg / 18.4 st (w/ Lithium-Ion battery)
Hill climbing ability: 17°
Maximum range on single charge: 45-60 km / 28-37.3 mi
Speed of response: Low
Maximum speed: 8 mph / 12 km/h
Acceleration system: Right-hand scooter throttle
Seat adjustment: Yes
Suspension: Full
Compliant with DVLA requirements: Yes
Scooter class: Class-3
Vehicle autonomy: Level 0
Optimal air humidity for storage and use: 30%
Working temperature of electronic parts: -10 °C to +40 °C
Turning radius: 220 cm / 86.6"

MOTOR

Power: 1000 W
Voltage: 60 V

BATTERY

Type: Lead-Acid Gel (can be upgraded to Lithium-Ion)
Capacity: 20 Ah
Voltage: 60 V
Life cycle: 600 full charges
Charge time: ≈8 h

BRAKES

Front: Hydraulic disc
Rear: Drum
Magnetic brake: No

WHEELS

Amount: 4
Type of tyres: Pneumatic tubeless
Size: F 3.00-8" / R 3.00-10"
Tyre pressure: 2.25 bar / 225 kPa / 32.5 PSI

LIGHTS

Low beam headlight: LED (integrated with lamp)
High beam headlight: LED (integrated with lamp)
Position lights: LED (integrated with lamp)
Tail lights: LED (integrated with lamp)
Front and rear signal lights: P21/5 W LED / R10W 12V
Brake light: P21/5 W LED
Hazard lights: P21/5 W LED / R10W 12V

FEATURES

Horn
Rearview mirrors
Steering wheel lock
Parking brake
Electric LED display
USB Port
Roof and windshield
Alarm
Remote control

FASTER +ROOF

